



Introducing Hawkin TruStrength, invented to provide you with invaluable data-confidence for empowered decision-making about your athletes and patients. This modular solution packs a big punch by providing precision biomechanical insights, allowing you to capture and view true data in real time. Hawkin TruStrength combines the aptitude of hand-held dynamometers and force plates but offers much more.

| | |
|----------------------|---|
| Size | 105 W X 80 H X 95 D mm |
| Weight | 1.7 lbs (800g) |
| Sensor Type | Load Cell |
| Resolution | 1 N |
| Capacity | 10,000 N / 1,000 kg Tensile and Compression |
| Sampling Rate | 1,200 Hz |
| Variables | Peak Force Mean Force Impluse Peak RFD Contact Duration |