




GEN 3

Hawkin TruStrength delivers data you can trust—anywhere, anytime. Designed for confident decision-making, this wireless modular system provides precise biomechanical insights in real time through the Hawkin App. Whether you're in the gym or the clinic, TruStrength empowers you with portable dynamometry—the difference is its fixed setup, which ensures the accuracy your athletes and patients deserve.

Size	92 x 79 x 67 mm (L x W x H)
Weight	1.95 lbs (0.88kg)
Sensor Type	S-Beam Strain Gauge
Precision	1 N
Range	+/-9,806N
Load Capacity	10,000 N
Sampling Rate	1,200 Hz
Battery Capacity	14.8Wh rechargeable lithium-ion (~150 Hours of battery life)
Metrics	Peak Force Average Force Average Net Force Total Impulse Net Impluse Peak RFD Duration Target Force Pretension